Quick Reference Guide

Easy•Touch®

Glucose Monitoring System

Display Messages

- HI** Meter detects blood glucose level higher than 600 mg/dL (33.3 mmol/L).
- LO** Meter detects blood glucose level lower than 20 mg/dL (1.1 mmol/L).
- HI** Ambient temperature was too high to perform the test.
- LO** Ambient temperature was too low to perform the test.
- Er 1** Test strip is not inserted properly/Defective test strip.
- Er 3** Meter malfunction. Please call customer service.
- L-0** Battery power is low. Meter will provide approximately 50 more measurements.
- Er 0** Battery power is low. Replace immediately.

Precautions for Use:

MHC Medical Products, LLC is not responsible for any injury resulting from incorrect use of the EasyTouch® Glucose Monitoring System and not in line with the instructions set forth in this operating manual.

Only use EasyTouch® Test Strip with the EasyTouch® Meter.

Any change or administer of medication based on the EasyTouch® blood glucose test results without the consent and advice of a physician or healthcare professional is not recommended.


BEFORE TESTING:

Read detailed instructions in your EasyTouch® User's Manual.

MHC Medical Products, LLC Questions or Comments?
8695 Seward Road, Call 877.358.4342
Fairfield, Ohio, 45011 www.EasyTouchProducts.com

To turn the Meter on, first insert the battery, then press the reset button located at the back of the Meter and proceed with the following instructions.

Test Summary


STEP 1

Insert a Test Strip, push all the way into test port until meter automatically turns on. A beep will sound and moving dotted lines will appear followed by a flashing "♦" symbol. Now you are ready to obtain a blood sample.

STEP 2

Place blood sample onto the end of the test strip. The meter will beep when enough blood has filled the strip, (0.5ul).

STEP 3

After 5 seconds, your blood glucose result will appear on the display and will automatically be stored into the Meter's memory.

Memory Feature

1. Press the **M** button to turn on the Meter. The date and time will appear first.
2. Press the **S** button once and the result of control solution test will be shown.*
3. Press the **S** button again to display your 7-day average, twice for 14-day average, and three times for 30-day average.*
4. Press the **S** button again and the meter will recall the last 300 test results in most recent order taken.*
5. When "000" appears, you have reached the end of your results.
6. Press the **M** button to turn off the meter.

*If there is no test results during this stage, the most recent test result will directly appear.

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MHC Medical Products, LLC Preguntas o comentarios?